



PicoPlus®

Highest Power & Versatility For Difficult Cases

- Highest Power of Any Multi-Wavelength Picosecond Laser
- Shorter Pulses (450ps) – Reduced Pain & Erythema
- Best of Nanosecond & Picosecond Pulse Durations
- Advanced Proprietary Technology & Elegant Ergonomics
- Proven Performance & Reliability - Hundreds Sold, Now Available in the USA!

Intelligently Designed with the Physician in Mind



INTUITIVE

Designed for the most busy user, the system is a dream to use: instantly switch between wavelengths and pulse durations without the need for recalibration.



INNOVATIVE

Industry exclusive 0.01 J/cm² settings permit safe and effective treatments on dark skin. The ultra-short picosecond pulses provide effective treatments without increased pain nor erythema.



EFFECTIVE

With high power for the most difficult cases, such as mild pigment (cafe-au-lait and freckles), recalcitrant tattoos, exotic ink colors, and more.



DEPENDABLE

Physicians around the world rely on PicoPlus for optimal outcomes. It self calibrates to maximize uptime.

Reproducible Results... Even on DIFFICULT Cases

WIDE RANGE OF APPLICATIONS

- Stubborn Pigment
- Post Inflammatory Hyperpigmentation
- Freckles
- Cafe-au-lait
- Nevus of Ota
- ABNOM
- Lentigos
- Seborrheic Keratoses
- Tattoos (Multicolor & Black)
- And More



Courtesy of G.S. Lee, MD



Courtesy of Dr. Hont, MD



Courtesy of B.S. Cho, MD

“The system offers many other benefits, including pigment removal, rejuvenation* and laser facial toning. The user-friendly interface comes with a toggle function for adjusting fluence energy settings very easily. Settings can be adjusted at any stage of treatment and the warm up time is very quick.”

Saleem Loghdey, MD

Exclusively Distributed in Benelux by:

LUTRONIC

Medi Invest Nv
Georges Wittouckstraat 20, 1600 Sint-Pieters-Leeuw, BELGIUM
www.medi-invest.eu - Info@medi-invest.eu - Tel:+32(0)2 381 2700

MEDI INVEST